ITENA

TOTAL C-RAM

AUTOMIX / DUAL-CURE

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Self-etch & self-adhesive resin cement for ceramic GB

PROPERTIES

Total C-RAM is a dual cured permanent adhesive resin cement specially formulated for ceramics cementation

INSTRUCTIONS FOR USE

1. Remove the provisional restoration. Note: Use only temporary cements that are eugenol - free. The use of eugenol-containing temporary luting materials can inhibit the polymerization process of Total C-RAM. 2. Rinse thoroughly. 3. Air dry the preparation (do not desiccate). 4. **Try in the restoration to ensure proper fit.** 5. Bleed cartridge before initial use. Place mixing tip on Auto-Mix syringe. 6. Dispensing the material can be done as follows: Inlay/Onlay, Mount the supplied intra oral tip onto the mixing tip and dispense the cement directly into the cavity preparation covering all surfaces. Seat the restoration gently onto the preparation allowing the cement to flow from all sides. Crown: Dispense the cement directly into the restoration or tooth preparation covering all surfaces. Seat the restoration gently onto the preparation allowing the cement to flow from all margins. Note: Material working time is 1– 2 minutes; Set time is 3 minutes (at oral temperature).. 7. When the restoration is properly seated, remove all excess cement. Excess cement is best removed in its gel state with a scaler or explorer. Gel state can be achieved by light curing excess with a light for 2 seconds or allowing to self-cure for 1.5-2 minutes after application. For best results, leave undisturbed for 1.5 minutes before light curing. 8. After removal of excess cement, light cure all surfaces including margins for 20 seconds. 9. Use finishing strips for proximal areas. Polish margins with polishing discs, cups, or points. 10. Check occlusion.

CAUTION

Uncured methacrylate resin may cause contact dermatitis and damage the pulp. Avoid contact with skin, eyes and soft tissue. Wash thoroughly with water after contact.
Note: Automix tips are designed for single patient use only.

STORAGE AND SHELF LIFE

Do not expose to temperature exceeding 24°C/ 75°F. Do not expose to direct sunlight. Refrigerate for maximum shelf life. Do not freeze. Shelf life is 24 months from date of manufacture.

RECOMMENDED CURE TIMES

L.E.Demetron I, 20 seconds or Optilux 501 in Boost mode 20 seconds, Ramp Mode 20 seconds or Regular Mode 20 seconds. For all other lights see manufacturer's recommendation.

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F : Conserver au sec
BG : Keep in a dry place

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